

TECHNICAL CHART

POLY ASP ULTRACLEAR EZ-EDGE AR

CYL (- ON -) -0.00 ~ -3.00					POWER	CYL (+ ON -) -0.00 ~ -2.00				
DIAMETER (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE		DIAMETER (mm)	E/T (mm)	C/T (mm)	BACK BASE	FRONT BASE
Φ72	2.0	2.0	3.85	3.83	0.00					
	2.0	1.7	4.07	3.83	0.25	3.83	3.65	2.0	1.8	Φ68
	2.3		4.30	3.83	0.50	3.83	3.43	2.2		
	2.6		4.52	3.83	0.75	3.88	3.20	2.3	1.6	
	2.9		4.75	3.83	1.00	3.88	2.97	2.5		
	3.2		4.98	3.83	1.25	3.88	2.75	2.6	1.4	
	3.5		5.20	3.83	1.50	3.88	2.53	2.8		
	3.4	3.61		1.75	3.88	2.33	2.6			
	3.6	3.38		2.00	3.88	2.10	2.9			
	3.9	3.16		2.25	4.20	2.22	2.9			
	4.2	1.4	2.93	2.50	4.20	2.00	3.1	1.0		
	4.5		3.16	2.75	4.20	1.78	3.4			
	4.8		2.93	3.00	4.20	1.56	3.6			
	5.0		5.65	2.70	3.25	4.20	1.35		3.9	
	5.2	2.48		3.50	5.28	2.20	4.0			
	5.5	2.26		3.75	5.28	1.98	4.3			
	5.8	2.03		4.00	5.28	1.76	4.5			
	5.9	1.2	6.10	1.80	4.25	5.55	1.80		4.0	1.0
	6.3			2.03	4.50	5.78		4.3		
	6.5			1.80	4.75	6.00	4.4			
6.8	1.58			5.00	5.78	1.35	4.6			
7.0	1.35			5.25	6.00		4.8			
7.3	1.13			5.50	6.22	0.90	5.0			
7.5	0.90			5.75	6.00		5.2			
7.7	0.68			6.00	6.22	5.4				
Φ72 CYL -0.00 ~ -2.00	8.3	1.2	6.56	0.90	6.25	1. "MINUS" LENS SPH AXIS: -0.00 ~ -6.00/-8.00 CYL AXIS: -0.00 ~ -3.00/-2.00 2. "PLUS" LENS SPH AXIS: +0.25 ~ +6.00 CYL AXIS: -0.00 ~ -2.00				
	8.5			0.68	6.50					
	8.8			0.46	6.75					
	9.0			0.23	7.00					
	9.5		7.01	0.46	7.25					
	9.7			0.23	7.50					
	10.2		7.47	0.46	7.75					
	10.4			0.23	8.00					