



SFSV 1.74 Tintable Primer Hard Coated Lens

				Tinting Rate Time (%)		
				10min	20min	30min
Best Guideline						
Surfacing	→ Backside Spin	H/C	→ Tint → AR Coating	40%	50%	65%
Surfacing	→ Tint	→ Backside Spin	H/C → AR Coating	25%	40%	52%

SFSV 1.74 TINTABLE LENS AVAILABILITY

BASE	DIAMETER (mm)	THICKNESS (mm)		RX RANGE
		CENTER	EDGE	
1.00	75	7.5	14.0	-16.00D to -6.25D
2.00		5.3	11.5	-6.00D to -2.75D
3.00		6.0	10.5	-2.50D to +1.00D
4.00		6.0	9.0	+1.25D to +2.00D
5.00		8.5	9.5	+2.25D to +3.50D
6.00		8.7	8.4	+3.75D to +5.00D
7.00		9.0	7.2	+5.25D to +6.50D
8.00		11.0	7.7	+6.75D to +8.00D